



Leman Instruments

High Technology Gas Generators

Three decades of experience in the field of high technology Gas Generators enable us to open Leman Instruments SAS, a company specialising in the creation, development, production and marketing of components, modules of scientific equipment involved in the analysis of solids, liquids and gases. We are focused in the production of standard products, OEMs, and in custom designed realisation. The company is located in Archamps Technopole, in the Geneva area.

Benefits:

- Reduces operating costs
- Return on investment in less than 2 years
- Improves resolution and detection limit of your instrument.
Provides high pressure stability
- Provides H₂, N₂ and Zero Air 24/7 at constant purity.
No contamination
- No more handling and storage of gas cylinders.
No more rental fees
- Lowest Footprint on the market
- Can be controlled remotely by PC or any supports
- Internal maintenance free

Control via touch panel or via iPad.

Technologies involved:

- Proton Exchange Membrane (PEM)
- Pressure Swing Absorption (PSA)
- Hot Catalysis (Oven).



SWISS TECHNOLOGIES & EUROPEAN MANUFACTURING



www.leman-instruments.com



Telephone: 0033 450394922 | Email: a.berteletti@leman-instruments.com